Colonoscopy versus FIT in Reducing Mortality from CRC (CONFIRM)

Impact of Mailing Time and Season on FIT Positivity



CONFIRM Study To	
	am
Study Co-Chairs.	
Jason A. Dominitz, MD, MPI	
West Lieven CCD Coordination Conter	
West Haven CSP Coordinating Center.	
Gary Jonnson, MS	
lassos Kyriakides, PhD	
Robert Wallace, ScD, MPH	
Albuquerque Pharmacy Coordinating Center:	
Katny Boardman, RPh	
Barbara Del Curto	
Other Executive Committee Members	
Dennis Ahnen, MD, Denver, CO	
Peter Guarino, PhD, Seattle, WA	
Tom Imperiale, MD, Indianapolis, IN	
David Lieberman, MD, Portland, OR	
Dawn Provenzale, MD, MS, Durham, NC	
Shahnaz Sultan, MD, Minneapolis, MN	
Aasma Shaukat, MD, Minneapolis, MN	







Current Status of Trial

- ~36 centers enrolling participants
- Recruitment milestones
- First randomization May 22, 2012
- 5,000th recruitment in June 2013
- 15,000th recruitment in April 2014
- 25,000th recruitment in March 2015
- 35,000th recruitment in February 2016

CONFIRM



Interventions: FIT

One sample

- OC-FIT CHEK analyzed with the Diana®
- Polymedco, Inc., Cortland Manor, NY
- Positivity: ≥20 µg hgb/g stool
- Centralized processing (Albuquerque, NM)
- Only samples arriving within 15 days of collection are processed

FIT Positivity by Season

If no collection date, specimen is processed

Samples are refrigerated upon arrival in lab

Mean Hgb Concentration by Season

Study	Spring	Summer	Fall	Winter	
Grazzini et al	27.6 ng/mL	25.2 ng/mL	29.2 ng/mL	29.5 ng/mL	<0.001
Cha et al [†]	0.45 ng/mL	0.25 ng/mL	0.20 ng/mL	0.32 ng/mL	<0.001
		Gut 2010; 59:1511-1515 Dig Dis Sci 2012; 57:2178-2183			

Study	Test	Spring	Summer	Fall	Winter	р
Grazzini et al	OC Sensor (20µg)	3.7%	3.5%	3.7%	4.2%	Not reported
Chausserie	OC Sensor (30µg)	2.9%	2.3%	3.0%	3.0%	0.159
Van Roon	OC Sensor (20 µg)		8.0%		9.7%	0.006
			Gut 20 Int J. C Am J G	10; 59:1511 ancer 2015 Sastroenterc	-1515 ;n136:1827-)/ 2012;107:	-1834 99-107

Transit Delay and FIT Positivity Population: Netherlands FIT: OC Sensor (10 μg hgb/g stool) Of the 6157 who returned a kit, 3767 reported the date of collection Transit Time N Positivity 14 days 3067 8.7% 25 days 705 6.0%

195

≥7 days

4.1%





Methods

- Initially, FIT mailed using first class mail
- Switched to US Priority Mail after discovery of some unacceptable shipping delays in some areas
- Chi-square test was used to compare positivity rate by season and the association between FIT positivity and shipping time





Results: Shipping Time (Days)

1 st Class (n=8020) 8.5 9.4 7	4-12	14	17		
Driarity Mail					
(N=16,141) 6.5 4.2 5	3-8	13	14		
p<0.0001					





	Association Between Season and FIT Positivity		
	Meteorological Season	Kits Processed (N)	Positivity (%)
j	Winter (Dec 1- Feb 28)	3754	7.3
	Spring (Mar 1- May 31)	5290	6.9
	Summer (June 1-Aug 31)	5712	5.8
	Fall (Sep 1- Nov 30)	5391	7.2
		p=0.007	



Summary

- Overall FIT positivity rate was 6.8%
- Positivity rates varied significantly by season, being lowest in the summer
- Shipping times of up to 14 days did not appear to affect FIT positivity rates
- The presence or absence of a sample collection date did not affect positivity rates

CONFIRM

Limitations

- We did not assess impact of variation in positivity rate on endoscopic findings, such as adenoma or advanced neoplasia detection
- Generalizability to FIT-based programs that are not nationally-based (with presumably smaller distances between participant and processing center) is are not clear



CONFIRM

Sample Return Time & Season with New Buffer

	Positivity
6130	4.1%
8643	4.1%
5598	4.6%
4797	4.3%
8488	4.1%
4818	4.5%
2268	3.9%
	6130 8643 5598 4797 8488 4818 2268

